Ultrasonic nebulizers and the enzymatic effect of rhDNase (Pulmozyme®), Extract (1)

Concluding report on in vitro investigations

Aim of investigations

Patients suffering from cystic fibrosis (CF) show an increase in the viscosity of the mucus produced in the lung due to a greatly elevated proportion of free DNA. This increase in viscosity leads to a deterioration in pulmonary function which can, however, be combated by various types of treatment.

Inhalation of rhDNase (Pulmozyme®) leads to a marked decrease in mucus viscosity and thus improved pulmonary clearance.

The aim of the investigations was to analyze enzymatic rhDNase activity before and after nebulization with the Multisonic® ultrasonic nebulizer manufactured by Schill, and to compare it to a compressed-air nebulizer (Pari-Master).

Materials and methods

During in vitro investigations, rhDNase (Pulmozyme®, manufactured by Roche) was nebulized using the Multisonic® ultrasonic nebulizer manufactured by Schill. Gel electrophoresis was used to determine enzyme activity before and 2.5 sec, 1 min, 5 min, 15 min, 30 min and 60 min after nebulization on high-molecular human DNA. The enzyme activity of non-nebulized enzyme and activity after nebulization with a compressed-air nebulizer (Pari-Master) before and 2.5 sec, 1 min, 5 min, 15 min, 30 min and 60 min after nebulization were also measured for comparison.

Discussion

In vitro investigations of the nebulizer types required for inhalational therapy with Pulmozyme® (Multisonic ultrasonic nebulizer and Pari-Master compressed-air nebulizer) have shown that there are no significant differences between the two of them as regards enzyme activity. Although different nebulization techniques are used, the effect on Pulmozyme® activity is the same for both systems and thus they both show a slight reduction in activity as compared to non-nebulized Pulmozyme®. In addition, the duration of nebulization does not appear to have an effect on enzyme efficacy for the two nebulizer types.

The extremely rapid fragmentation of long-chain DNA and thorough breaking up into pieces after 60 minutes is striking in the case of the non-nebulized enzyme. In these areas, Pulmozyme® efficacy of both nebulizer types is slightly weaker than that of the non-nebulized enzyme. Heat developing during operation of the ultrasonic nebulizer could be responsible for this, as could a sonification effect caused by ultrasound.


On the other hand, shearing forces occurring in the Pari-Master compressed-air nebulizer could slightly affect enzyme activity. A change in reactivity could also occur due to the nebulization process and the close contact of the enzyme with air or O₂ which occurs in both devices.

Despite this, the efficacy of Pulmozyme® is only slightly reduced by the use of the two nebulizers tested. Other factors such as the uniform distribution of the aerosol in the lung could be more significant in this context and should be considered in future investigations.

Summary


In the in vitro investigations conducted on the enzymatic efficacy of Pulmozyme® after short and long nebulization periods with the Multisonic ultrasonic nebulizer, no significant difference could be found as compared to nebulization with the Pari-Master compressed-air nebulizer. A slight reduction in efficacy of the nebulized enzyme as compared to the non-nebulized enzyme may be observed but only, however, for certain DNA lengths.

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